ADLITETM

- Sealed concentric beam, disposable diagnostic penlight
- Activated by depressing pocket clip
- Available with or without pupil gauge
- Measures 5" long, 1/2" diameter
- Sold in 6 pack



ADLITETM

AMERICAN DIAGNOSTIC CORPORATION

Diagnostic Penlights

Ordering Information

ADLITE™ Disposable Penlight, 6/pk

#351P ADLITE™ Disposable Penlight with Pupil Gauge, 6/pk

ADLITETM PLUS

- Sealed concentric beam, disposable diagnostic penlight
- Activated by depressing pocket clip
- · White available with or without pupil gauge
- Measures 5" long, 1/2" diameter
- Sold individually or in display cylinder of 22
- Available in 7 colors



Ordering Information

#356* ADLITE™ Disp. Penlight, specify color

#358TVB Disp. Penlight, 22/Cylinder - Blk, Purple, Teal





55 Commerce Drive Hauppauge, New York 11788 tel: 1-800-232-2670 fax: 631-273-9659

www.adctoday.com

REUSABLE PENLIGHTS

Three distinctly different models to suit your personal preferences.



ADLITE PROTM

Our finest pocket light, the ADLITE™ Pro is the physician's choice for evaluation of the pupilary reflex.

- LED lamp for bright, white illumination for true tissue rendering.
- LED lamp provides 10,000 hour use eliminating the high cost of replacement lamps.
- Durable CNC machined brass housing in a striking black matt finish.
- Positive feel power switch
- Measures 5%" long, 1/2" diameter
- · Convenient Pocket Clip
- Complete with 2 "AAA" batteries and an attractive gift box.
- · Lifetime warranty on lamp cartridge, 2 year warranty on penlight.

ORDERING INFORMATION

#355BK Adlite Pro™, Physician's Penlight



ADLITE IITM

- · Plastic housing for lightweight comfort
- Pocket clip activated
- Measures 5" long by 1/2" diameter
- Two "AAA" batteries included

ORDERING INFORMATION

DIAGNOSTIC PENLIGHT

ADLITE II™ Reusable Penlight

METALITETM

- Brushed aluminum housing
- Locking push toggle activated
- Measures 5 ¼" long by 1/2" diameter
- Two "AAA" batteries included

ORDERING INFORMATION

METALITE™ Reusable Penlight



Display packaging available. Add "Q" suffix to order number.